

AAO 2018 ART + SCIENCE

Exhibitor Guide and Hall Map

Presented by EyeNet® Magazine



Protecting Sight. Empowering Lives.®

BUILT FOR THE NEXT GENERATION. BUT READY FOR THIS ONE.



0

) (Q)

0

0

0

BAUSCH+LO

Experience Stellaris Elite[™]–where real-time responsiveness enables exceptional control.

- V ----





AAO Booth 3044 (800) 338-2020 | StellarisElite.com

Indications: The Bausch + Lomb Stellaris Elite[®] vision enhancement system is intended for the emulsification and removal of cataracts, anterior and posterior segment vitrectomy. The system is designed for use in both anterior and posterior segment surgeries. It provides capabilities for phacoemulsification, coaxial and bimanual irrigation/aspiration, bipolar coagulation, vitrectomy, viscous fluid injection/removal and air/fluid exchange operations. The Stellaris Elite vision enhancement system configured with the laser module is additionally intended for retinal photocoagulation and laser trabeculoplasty. Contraindications: All Systems: Use of accessories not designated by Bausch + Lomb for use with this equipment may result in serious permanent patient injury, adverse surgical outcome, or damage to the equipment. Systems with Laser Module: Photocoagulation is not indicated for patients without pigmentation (albino eyes). In addition, Laser Indirect Ophthalmoscope (LIO) is not indicated for cases involving laser photocoagulation within the arcades. Warnings: All Systems: Implantable defibrillators present a risk of injury if triggered by a fibrillatory event during intraocular surgery. Electromagnetic interaction between the phacoemulsification (phaco) handpiece and an implanted cardiac pacemaker is unlikely, but cannot be ruled out. Systems with Laser Module: All support personnel who are present during laser treatment must wear appropriate laser protective eyewear. **DO NOT** look directly into the aiming or treatment laser beam. Use of unapproved delivery devices may cause inaccurate laser delivery which could result in serious permanent patient injury. When using the Vitesse[®] mandpiece: Use only the Entry Site Alignment (ESA) devices provided with the Vitesse of the laser system may lead to dangerous situations and severe injuries. All Systems: See the appropriate Operator Manual for detailed directions, proper use, and full risk and safety information. See Vitesse handpiece pack Instructions for Uses fon

[®]/[™] are trademarks of Bausch & Lomb Incorporated or its affiliates. All other brand/product names are trademarks of their respective owners. ©2017 Bausch & Lomb Incorporated. SEP.0051.USA.17

Aerie Pharmaceuticals, Inc.

Booth 747

Ellex

Aerie is an ophthalmic pharmaceutical company focused on the discovery, development and commercialization of first-in-class therapies for the treatment of patients with open-angle glaucoma, retina diseases and other diseases of the eye. We have a passion for finding novel ways to protect vision – beginning with the first new drug class introduced for glaucoma in 22 years. Learn more at www. aeriepharma.com.



Visit Ellex to see the latest in laser and ultrasound technologies for the diagnosis and treatment of cataract, glaucoma, and retinal conditions, including Laser Floater Treatment (LFT). Our solutions now include ABiC with iTrack™ for minimally invasive glaucoma surgery (MIGS). We will also have exhibit presentations on Laser Floater Treatment, ABiC with the iTrack microcatheter, and

Retinal Rejuvenation with 2RT.



61-8-7675-784 www.ellex.com

Ocular Instruments, Inc.

Booth 2619

Booth 1044



We invite you to see our latest innovations for MIGS, Laser, Surgical and Diagnostic procedures. Featuring the NEW Mainster Wide Field EX and SecureFlex® Surgical Gonio and lenses.

For more information on our complete line of superior quality products visit our website, better yet stop by Booth 2619 and take advantage of our 30 day no risk trial!

425-455-5200 www.ocularinc.com

OMIC - Ophthalmic Mutual Insurance Company Booth 1004



Visit the OMIC Exhibit and see our new sample patient education materials, informed consent documents, and other loss prevention tools.

OMIC is the largest insurer of ophthalmologists in the United States and is endorsed by more than 50 ophthalmic state and subspecialty societies, including the exclusive sponsorship of the American Academy of Ophthalmology.

800-562-6642 www.omic.com

855-237-4379 www.aeriepharma.com

Bausch + Lomb

The Stellaris Elite[™] vision enhancement system is a next-generation surgical platform from Bausch + Lomb. With a comprehensive combination of leading edge technologies, Stellaris Elite is highly responsive at every moment of your procedure to enable exceptional stability. And only Stellaris Elite offers the vitreous removal capabilities of the Vitesse[™] hypersonic vitrectomy system and the chamber stabilization system. To reserve your immersive wet lab experience with Bausch + Lomb's new technologies, please visit booth #3044.



Booth 3933

Booth 3044

800-553-5340 www.bausch.com

Corne<u>aGen</u>



Geuder Glass Cannula for DMEK

CorneaGen is a mission-driven company committed to eliminating corneal blindness and to transforming how corneal surgeons treat and care for the cornea. We have a growing family of innovative products including the **Pre-Loaded Geuder Glass Cannula for DMEK** and the **KAMRA* Inlay** corneal implant. Visit us to learn more about how we are innovating the next generation of cornea care.

206-701-5840 www.corneagen.com/chicago



elcome to Chicago! *EyeNet's Exhibitor Guide* is designed to maximize your time on the show floor. The guide's alphabetical business listings help you quickly identify companies to visit and highlight some of their new products and services. Use the foldout map to navigate the exhibit floor, and use the notes page in the back to plot your course through the hall or to take notes during vendor meetings. We hope your time at AAO 2018 is productive and enjoyable.

Sincerely

Dale E. Fajardo, EdD, MBA Publisher, *EyeNet Magazine*



AMERICAN ACADEMY OF OPHTHALMOLOGY® Protecting Sight. Empowering Lives. **DISCLAIMER:** The product information and claims provided in the profiles included in *EyeNet's Guide to Academy Exhibitors* are those of the manufacturers and have not been verified, nor does the appear-

ance of a product constitute an endorsement by EyeNet or the American Academy of Ophthalmology. Information in this guide is current as of Aug. 22, 2018. Please check aao.org for updates to Exhibit Hall booths. American Academy of Ophthalmic Executives[®], *EyeNet*[®], EyeSmart[®], Foundation of the American Academy of Ophthalmology[®],

IRIS[®] Registry, ISRS[®], Museum of Vision[®], ONE[®], and Preferred Practice Patterns[®], are trademarks of the American Academy of Ophthalmology[®]. All other trademarks are the property of their respective owners.



66 Vision Tech Co., Ltd	Booth 830	Advice Media Booth 765
86-5126-7770-255 www.66vision.com		435-575-7470 www.advicemedia.com
A Fashion Hayvin, Inc.	Booth 1803	Aeon Astron Europe B.V. Booth 4526
321-239-0184 www.conventionjewelry.com		31-7133-2-2280 www.aeonastron.com
A.R.C. Laser GmbH 49-91121-779-0 www.arclaser.com	Booth 2453	Aerie Pharmaceuticals, Inc. Booth 747
AA Vision, Inc. 91-1126-41-5788 www.aavisioninc.com	Booth 1603	Aerie is an ophthalmic pharma- ceutical company focused on the discovery, development and commercialization of first-in-class therapies for the treatment of patients with open-angle glau- coma, retina diseases and other diseases of the eye. We have a
AAO EyePlay Experience 415-561-8500 www.aao.org	Booth 2581	passion for finding novel ways to protect vision - beginning with the first new drug class intro- duced for glaucoma in 22 years.
Aaren Scientific, Inc. 909- 906-5400 www.aarenlabs.com	Booth 2267	Learn more at www.aeriepharma. com. 855-237-4379 www.aeriepharma.com
AbbVie 800-255-5162 www.abbvie.com	Booth 1049	AJL OphthalmicBooth 325334 945-298-256 www.ajlsa.com
Accutome 800-979-2020 www.accutome.com	Booth 2244	Akorn, Inc.Booth 3700800-932-5676 www.akorn.com
Action Medical Pvt Ltd 91-2137-25-3617 www.actionmedihealth.com	Booth 416	AKtive SrlBooth 294839-06-8864-0509 www.aktive-corp.com
Acuity Pro, Inc. 580-243-1301 www.acuitypro.com	Booth 3167	AL.CHI.MI.A. Srl Booth 2512 39-0498-96-2074 www.alchimiasrl.com
Addition Technology, Inc. 34-945-298-256 www.ajlsa.com	Booth 3253	Alcon Booth 2708 817-293-0450 www.alcon.com
Advanced Data Systems 800-899-4237 www.adsc.com	Booth 659	Alex Surgical Co., LtdBooth 357486-2585-9967-18 www.alexsurgical.com
AdvancedMD 801-984-9500 www.advancedmd.com	Booth 4605	Alimera Sciences Booth 1535 678-990-5740 www.alimerasciences.com

Allergan, Inc. Booth 1308 714-246-4500 www.allergan.com	American Academy of Ophthalmology (AAO) Booth 508
	AMERICAN ACADEMY OF OPHTHALMOLOGY*
Allergan Medical Affairs Booth 1800	Protecting Sight. Empowering Lives.
714-246-4500 www.allergan.com	Find Resources to Enhance Your Clinical and Business Success at the Academy Resource Center
Alphaeon Credit Booth 1947	 Save 10% on products. Demo and purchase clinical education, patient education, and practice management products.
920-360-1794 www.myalphaeoncredit.com	 Discover how the Academy's ONE® Network and IRIS® Registry can improve the way you practice. Attend Special Events—see the complete listing at aao.org/
ALT Medical Device Booth 3164	resourcecenter.
786-657-3771 www.mdalt.com	415-561-8500 www.aao.org
Altaire Pharmaceuticals, Inc. Booth 2101	American Express Booth 238
631-722-5988 www.nanotearseyedrops.com	212-640-2000 www.open.com
AMA Optics, Inc. Booth 250	American Society of Cataract Booth 2035 & RefractiveSurgery (ASCRS)
877-744-3937 www.amaoptics.com	703-591-2220 www.ascrs.org
Ambler Surgical Booth 2148	American Society of Ophthalmic Booth 2035 Administrators (ASOA)
610-280-7361 www.amblersurgical.com	703-591-2220 www.asoa.org
AMCON Booth 3304	Amring Pharmaceuticals, Inc. Booth 2278
314-849-7773 www.crestpointmgt.com	610-804-4663 www.amringusa.com
American Academy of Booth 508 Ophthalmic Executives (AAOE)	Angiotech Booth 2845
	877-991-1110 www.surgicalspecialties.com
AMERICAN ACADEMY OF OPHTHALMOLOGY* Protecting Sight. Empowering Lives.	Anodyne Surgical Booth 2138
Tackle Your Business Challenges With Expert Advice	636-447-1010 www.anodynesurgical.com
Visit the Academy Resource Center for help from AAOE experts: • Schedule a free, one-on-one, 20-minute consultation with a	AO Lab - American Ophthalmic Lab Booth 1600
 practice management specialist. Gain a better understanding of MACRA and MIPS. Find the answers to your coding conundrums. 	410-647-5140 www.aolab.com
 Get the newest instructional and reference materials to help you improve your skills and increase profitability. 	CHECK YOUR EMAIL FOR AAO 2018 DAILY
415-561-8500 aao.org/practice-management	For concise summaries of each day's clinical highlights, read AAO 2018 Daily, a brief bulletin emailed in the evening to you and to more
USE THE MOBILE MEETING GUIDE	than 70,000 of your colleagues in the global ophthalmic comu-

nity. AAO 2018 Daily is reported onsite in Chicago Oct. 26-29, and

of news from Subspecialty Day and AAO 2018. Additionally, these

articles will be posted online at aao.org/eyenet/daily.

allows ophthalmologists at the meeting and at home to stay on top

USE THE MOBILE MEETING GUIDE

The Mobile Meeting Guide is available at aao.org/mobile and contains all meeting and course information including abstracts, posters, and videos.

Mastering intraocular pressure control with Rhopressa®



What makes once-daily Rhopressa® different

ÛÛ	Consistent mean IOP reduction of 5 mmHg in patients across a range of baseline IOPs ¹	- <u>``</u>	Once-daily dosing to simplify dosing regimens ¹
0	Mild ocular adverse events and no known contraindications opens up treatment options ¹	Ó	Unique mechanism of action for patients who may benefit from improved trabecular aqueous outflow ¹

Visit Rhopressa.com to learn more about this innovative IOP-lowering treatment.

Share the Vision with Aerie Pharmaceuticals—join us at AAO Chicago, Booth 747

INDICATION AND IMPORTANT SAFETY INFORMATION

INDICATION

Rhopressa® (netarsudil ophthalmic solution) 0.02% is a Rho kinase inhibitor indicated for the reduction of elevated intraocular pressure in patients with open-angle glaucoma or ocular hypertension.

Dosage and Administration: The recommended dosage is one drop in the affected eye(s) once daily in the evening.

IMPORTANT SAFETY INFORMATION

Dosage and Administration: Twice a day dosing is not well tolerated and is not recommended. If Rhopressa® is to be used concomitantly with other topical ophthalmic drug products to lower IOP, administer each drug product at least 5 minutes apart.

Warnings and Precautions:

Bacterial Keratitis: There have been reports of bacterial keratitis associated with the use of multiple-dose containers of topical ophthalmic products. These containers had been inadvertently contaminated by patients who, in most cases, had a concurrent corneal disease or a disruption of the ocular epithelial surface.

Adverse reactions: The most common ocular adverse reaction observed in controlled clinical studies with Rhopressa® dosed once daily was conjunctival hyperemia which was reported in 53% of patients. Other common (approximately 20%) adverse reactions were: corneal verticillata, instillation site pain, and conjunctival hemorrhage. Instillation site erythema, corneal staining, blurred vision, increased lacrimation, erythema of eyelid, and reduced visual acuity were reported in 5-10% of patients. The corneal verticillata seen in Rhopressa®-treated patients were first noted at 4 weeks of daily dosing. This reaction did not result in any apparent visual functional changes in patients. Most corneal verticillata resolved upon discontinuation of treatment.

Please see brief summary of full Prescribing Information on the adjacent page.

Reference:

1. Rhopressa Prescribing Information. Irvine, CA: Aerie Pharmaceuticals, Inc; 2017.



www.aeriepharma.com

©2018 Aerie Pharmaceuticals, Inc. All rights reserved. Durham, NC 27703

RHOPRESSA® (netarsudil ophthalmic solution) 0.02% Rx Only

BRIEF SUMMARY

Consult the Full Prescribing Information for complete product information.

INDICATIONS AND USAGE

RHOPRESSA* (netarsudil ophthalmic solution) 0.02% is indicated for the reduction of elevated intraocular pressure (IOP) in patients with open-angle glaucoma or ocular hypertension.

DOSAGE AND ADMINISTRATION

The recommended dosage is one drop in the affected eye(s) once daily in the evening.

If one dose is missed, treatment should continue with the next dose in the evening. Twice a day dosing is not well tolerated and is not recommended. If RHOPRESSA is to be used concomitantly with other topical ophthalmic drug products to lower IOP, administer each drug product at least 5 minutes apart.

WARNINGS AND PRECAUTIONS

Bacterial Keratitis

There have been reports of bacterial keratitis associated with the use of multiple-dose containers of topical ophthalmic products. These containers had been previously contaminated by patients who, in most cases, had a concurrent corneal disease or a disruption of the ocular epithelial surface.

Use with Contact Lenses

RHOPRESSA contains benzalkonium chloride, which may be absorbed by soft contact lenses. Contact lenses should be removed prior to instillation of RHOPRESSA and may be reinserted 15 minutes following its administration.

ADVERSE REACTIONS

Clinical Trials Experience

Because clinical studies are conducted under widely varying conditions, adverse reaction rates observed in the clinical studies of a drug cannot be directly compared to rates in the clinical studies of another drug and may not reflect the rates observed in practice.

The most common ocular adverse reaction observed in controlled clinical studies with RHOPRESSA dosed once daily was conjunctival hyperemia which was reported in 53% of patients. Other common (approximately 20%) ocular adverse reactions reported were: corneal verticillata, instillation site pain, and conjunctival hemorrhage. Instillation site erythema, corneal staining, blurred vision, increased lacrimation, erythema of eyelid, and reduced visual acuity were reported in 5-10% of patients.

Corneal Verticillata

Corneal verticillata occurred in approximately 20% of the patients in controlled clinical studies. The corneal verticillata seen in RHOPRESSA-treated patients were first noted at 4 weeks of daily dosing. This reaction did not result in any apparent visual functional changes in patients. Most corneal verticillata resolved upon discontinuation of treatment.

USE IN SPECIFIC POPULATIONS

Pregnancy

There are no available data on RHOPRESSA use in pregnant women to inform any drug associated risk; however, systemic exposure to netarsudil from ocular administration is low. Intravenous administration of netarsudil to pregnant rats and rabbits during organogenesis did not produce adverse embryofetal effects at clinically relevant systemic exposures.

Animal Data

Netarsudil administered daily by intravenous injection to rats during organogenesis caused abortions and embryofetal lethality at doses $\geq 0.3 \text{ mg/kg/day}$ (126-fold the plasma exposure at the recommended human ophthalmic dose [RHOD], based on C_{max}). The no-observed-adverse-effect-level (NOAEL) for embryofetal development toxicity was 0.1 mg/kg/day (40-fold the plasma exposure at the RHOD, based on C_{max}).

Netarsudil administered daily by intravenous injection to rabbits during organogenesis caused embryofetal lethality and decreased fetal weight at 5 mg/kg/day (1480-fold the plasma exposure at the RHOD, based on C_{max}). Malformations were observed at \geq 3 mg/kg/day (1330-fold the plasma exposure at the RHOD, based on C_{max}), including thoracogastroschisis, umbilical hernia and absent intermediate lung lobe. The NOAEL for embryofetal development toxicity was 0.5 mg/kg/day (214-fold the plasma exposure at the RHOD, based on C_{max}).

Lactation

There are no data on the presence of RHOPRESSA in human milk, the effects on the breastfed infant, or the effects on milk production. However, systemic exposure to netarsudil following topical ocular administration is low, and it is not known whether measurable levels of netarsudil would be present in maternal milk following topical ocular administration. The developmental and health benefits of breastfeeding should be considered along with the mother's clinical need for RHOPRESSA and any potential adverse effects on the breastfed child from RHOPRESSA.

Pediatric Use

Safety and effectiveness in pediatric patients below the age of 18 years have not been established.

Geriatric Use

No overall differences in safety or effectiveness have been observed between elderly and other adult patients.

NONCLINICAL TOXICOLOGY

Carcinogenesis, Mutagenesis, Impairment of Fertility

Long-term studies in animals have not been performed to evaluate the carcinogenic potential of netarsudil. Netarsudil was not mutagenic in the Ames test, in the mouse lymphoma test, or in the *in vivo* rat micronucleus test. Studies to evaluate the effects of netarsudil on male or female fertility in animals have not been performed.



Manufactured for: Aerie Pharmaceuticals, Inc., Irvine, CA 92614, U.S.A.


For more information, go to www.RHOPRESSA.com or call 1-855-AerieRx (1-855-237-4379).

RHOPRESSA is a registered trademark of Aerie Pharmaceuticals, Inc.

U.S. Patent Nos.: 8,450,344; 8,394,826; 9,096,569; 9,415,043

Appasamy Associates Booth 744	Aurora Surgical, LLC Booth 2611
40 years of dedicated manufacturing of ophthalmology products. 25 years of participation at AAO.	727-821-3303 www.aurorasurgical.com
Ophthalmic Refraction Units; Slitlamps; Applanation Tonometer; Slitlamp Digital Imaging System; Auto Perimeter; Manual Refrac- tor; Auto Ref-Keratometer; LED HD Visual Acuity Charts; Indirect Ophthalmoscopes; Ultrasound A-Scan/B-Scan/UBM/Pachymeter; Ophthalmic Lasers; Operating Microscopes LED; Phaco/Vitrecto-	Avedro, Inc.Booth 720844-528-3376 www.avedro.com
my Systems; Mobile Operating Table; Hydrophilic IOLs; Preloaded Hydrophobic IOLs; Ophthalmic Viscoelastic Devices; Micro Surgi- cal Instruments/Knives.	Avellino Labs Booth 437 650-396-3741 www.avellinolab.com/us
91-44-30101423 www.appasamy.com	AZTACH Medical Device, Inc. Booth 2678 909-755-5533 www.aztach.com
ArbiMed Inc.Booth 2477832-461-1752 www.arbimed.com	Bankers Healthcare GroupBooth 2874866-297-4664 www.bankershealthgroup.com
ARCADOPHTABooth 275633-5-614-0523-5 www.arcadophta.com	Bausch + LombBooth 3044The Stellaris Elite™ vision enhancement system is a next-generation surgical platform from Bausch + Lomb. With a comprehensive
ArcScan, Inc.Booth 4463877-363-7226 www.arcscan.com	combination of leading edge technologies, Stellaris Elite is highly responsive at every moment of your procedure to enable exceptional stability. And only Stellaris Elite offers the vitreous removal capabilities of the
Army Medical Recruiting Booth 1700 888-550-2769 www.goarmy.com	Vitesse [™] hypersonic vitrectomy system and the chamber stability of Adaptive Fluidics [™] chamber stabilization system. To reserve your immersive wet lab experience with Bausch + Lomb's new technologies, please visit booth #3044.
Asia-Pacific Academy of OphthalmologyBooth 12028523943-5827 www.apaophth.org	800-553-5340 www.bausch.com
ASICO, LLC Booth 2604 630-986-8032 www.asico.com	Beauty and BeyondBooth 968, 2764702-883-0945 www.truffoire.com
Association for Research in Vision and Ophthalmology (ARVO)Booth 2658240-221-2900 www.arvo.org240-221-2900 www.arvo.org	Beaver - Visitec InternationalBooth 3356866-906-8080 www.bvimedical.com
AST Products, Inc. Booth 1949 978-667-4500 www.astp.com	Beneficial Equipment Finance Corp.Booth 3371877-880-9020 www.thebeneficial.com
AurolabBooth 194491-452-3096-100 www.aurolab.com	Benjamin Biomedical, Inc.Booth 3573727-343-5503 www.benjaminbiomedical.com

Besse MedicalBooth 340800-543-2111 www.besse.com	Canon U.S.A., Inc.Booth 3422631-330-5000 www.usa.canon.com/eye-care
Beyeonics Surgical, LtdBooth 1281972-48-318-899 www.elbitsystems.com	CareCloudBooth 1270305-265-4200 www.carecloud.com
BionikoBooth 648507-246-6456 www.bioniko.com	CareCreditBooth 2728800-300-3046 www.carecredit.com
Bio-Tissue, Inc.Booth 2736888-296-8858 www.biotissue.com	Cassini TechnologiesBooth 2670888-660-6965 www.i-optics.com
BlephEX, LLCBooth 337248-880-8753 www.rysurg.com	Ceatus Media GroupBooth 1845858-454-5505 www.ceatus.com
Blythe MedicalBooth 2573706-654-3210 www.blythemedical.com	Centervue, Inc.Booth 645408-988-8404 www.centervue.com
BreviumBooth 1702801-854-5400 www.brevium.com	CheckedUpBooth 652305-722-5498 www.checkedup.com
Bruder Healthcare CompanyBooth 1903888-827-8337 www.bruder.com	Chona Surgical, Co.Booth 235991-9582-21-5800 www.chonasurgicalco.com
Bryn Mawr Communications, LLCBooth 2803484-581-1800 www.bmctoday.com	Chongqing Kanghua Ruiming S & T Co., Ltd Booth 965 86-2363-2273-33 www.kanghuaa.com
BSM ConsultingBooth 350775-832-0600 www.bsmconsulting.com	Chongqing Sunkingdom MedicalBooth 2470Instruments Co., Ltd86-2368-6439-90 www.skoptic.net
C.S.O. Srl Booth 2935 39-0557-2219-40 www.csoitalia.it	Cilita, Ltd Booth 2411 7-49127-796-78 www.cilita.com
CalTray Booth 2403 413-776-1328 www.caltray.com	CIMA Technology Booth 3222 724-733-2627 www.cimatechinc.com
Cambrian Medical, Inc.Booth 2256408-996-0205 www.cambrianmed.com	Clearside Biomedical, Inc.Booth 3274678-270-3631 www.clearsidebio.com





Designer Drugs PharmacyBooth 451423-954-2585 www.designerdrugs.us	Duckworth & Kent LtdBooth 3900314-849-7773 www.crestpointmgt.com
Designs For Vision, Inc.Booth 2900800-345-4009 www.designsforvision.com	Duke Eye Center Booth 762
Device OpticalBooth 2865305-888-9499 www.deviceoptical.com	Duke Eye Center
Dexta CorporationBooth 2150707-255-2454 www.dexta.com	Duke Eye Center is one of the leading eye institutions in the world, and for almost 30 years, has been consistently ranked in the Top 10 by the US News and World Report Best Hospitals rankings. Our internationally-renowned clinical and research faculty are dedicated to curing eye disease worldwide though excellence in research, education, and
DGH Technology, Inc. Booth 2204 610-594-9100 www.dghtechnology.com	patient care. 919-668-1345 www.dukeeyecenter.duke.edu
Diagnosys LLCBooth 3021978-458-1600 www.diagnosysllc.com	Dutch Ophthalmic, USABooth 1335800-753-8824 www.dorc.eu
Diamatrix, LtdBooth 2104281-367-8081 www.diamatrix.com	Eagle LabsBooth 1940909-481-0011 www.eaglelabs.com
DigiSight Technologies, Inc.Booth 953650-223-5560 www.digisight.net	Eckert WordellBooth 2903269-388-7313 www.eckert-wordell.com
Digital Heat CorporationBooth 3576800-970-2904 www.heatedeyepad.com	Eidolon Optical, LLCBooth 730781-400-0586 www.slitlamp.com
Diopsys, Inc.Booth 1547973-244-0622 www.diopsys.com	Einstein MedicalBooth 2610866-308-8553 www.einsteinmedical.com
Doctorsoft CorporationBooth 251855-377-6381 www.doctorsoft.com	READ AAO 2018 NEWS Look for <i>AAO 2018 News</i> in 3 locations: near the escalator on Level 1; next to the bridge in Lakeside; and in the Resource Center, Booth 508.
Dongguan HEC Pharmaceutical R&D Co., Ltd Booth 2773 86-7698-5315-888 www.hecpharm.com/en/index.aspx	Subspecialty Day edition. The Subspecialty Day program directors recommend presentations to attend. Learn about ergonomics in your office. In addition, browse the year's most-read "Journal Highlights."
DORC International, bvBooth 1335800-753-8824 www.dorc.eu	AAO 2018 edition. Pick up the second edition for an in-depth profile of the 2018 Laureate, Steven T. Charles, MD. Other highlights include the Academy President's Guests of Honor, and a feature on the Museum of Vision's exhibit on early ophthalmic medicine.

Ellex Booth 1044 Visit Ellex to see the latest in laser and ultrasound technologies for the diagnosis and treatment of cataract, glaucoma, and retinal conditions, including Laser Floater Treatment (LFT). Our solutions now include ABiC with iTrack™ for minimally	European Society of Ophthalmology (SOE)Booth 290244-78302-2-1032 www.soevision.orgEuropean Society of Retina SpecialistsBooth 1405(EURETINA)353-12192 www.euretina.org
with iTrack [™] for minimally invasive glaucoma surgery (MIGS). We will also have exhibit presentations on Laser Floater Treatment, ABiC with the iTrack microcatheter, and Retinal Rejuvenation with 2RT. 800-824-7444 www.ellex.com	Excellent Hi-Care Pvt LtdBooth 25291-4423-6303-85 www.excellenthicare.comExhibitor LocatorExhibitor LocatorBooth 3400415-561-8500 www.aao.orgExhibitor Locator
Ellman, A Cynosure Company Booth 2253 516-594-3333 www.ellman.com	Eye Care and Cure Booth 2270 800-486-6169 www.eyecareandcure.com Double 1200
ELSEVIERBooth 3200215-239-3900 www.elsevier.comEnhanced Medical ServicesBooth 4257888-781-4458 www.emseye.com	Eye Care LeadersBooth 4308855-691-0126 www.eyecareleaders.comEye Designs, LLCBooth 3453800-346-8890 www.eyedesigns.com
Enhanced VisionBooth 4314800-440-9476 www.enhancedvision.com	Eye Photo Systems, Inc.Booth 4305503-488-5476 www.eyephotosystems.com
Epsilon EyeCare, USABooth 2400909-984-4614 www.epsilonusa.comEschenbach Optik	EyeKon MedicalBooth 2852727-793-0170 www.eyekonmedical.comEyekor, Inc.Booth 1550
800-487-5389 www.eschenbach.com Essilor Instruments Booth 1134 855-393-4647 www.essilorinstrumentsusa.com	608-206-3614 www.eyekor.comEyeMD EMR Healthcare Systems, Inc.Booth 756877-239-3367 www.eyemdemr.com
ESW – Vision Booth 826 33-1-348-7064-9 www.esw-vision.com	Eyenovia Booth 4266 917-289-1117
European Society of Cataract & RefractiveBooth 1403Surgeons (ESCRS)353-120-9-1100 www.escrs.org	Eyeol UK LimitedBooth 294644-15824-7700-2 www.eyeol.co.uk

EyeOn SwissFlexBooth 2873954-284-5663 www.swissflex-usa.com	Focus ROP, LLCBooth 4459610-955-5995 www.focusrop.com
Eyepoint Pharmaceuticals Booth 2775 617-926-5000 www.eyepointpharma.com	Focus Vision SupplementsBooth 3072703-832-5869 www.focusvitamins.com
EyeSys Vision, Inc. Booth 2511 281-885-3800 www.eyesys.com	Freedom Ophthalmic Pvt LtdBooth 31691-9585-5993-01 www.freedomophthalmic.com
Eyevance Pharmaceuticals LLC Booth 2869 919-724-7577 www.eyevance.com	Fresnel Prism and Lens Co.Booth 1604800-544-4760 www.fresnel-prism.com
EyeWorld Booth 1835 703-591-2220 www.eyeworld.org	GENENTECH Booth 4026 650-225-1000 www.gene.com
Fashion Optical Displays Booth 4523 800-824-4106 www.fashionoptical.com	Geuder AG Booth 3900 314-849-7773 www.crestpointmgt.com
FB Eye Equipments Booth 137 305-787-8366 www.fbeyeequipments.com	Glacial Multimedia Booth 3063 207-878-5900 www.glacial.com
FCI Ophthalmics Booth 3435 800-932-4202 www.fci-ophthalmics.com	Glaucoma Research FoundationBooth 1677415-986-3162 www.glaucoma.org
Feather Safety Razor Co., LtdBooth 250281-664-58-1638 www.feather.co.jp	Glaukos Booth 1345 949-367-9600 www.glaukos.com
First Insight Corporation Booth 2103 2015 ONC Certified. Easily customize exam templates to fit	Global Ophthalmic Pvt LtdBooth 116891-9444-66-6400 www.globalopthalmic.com
workflows. MIPS alerts and calculators, AAO IRIS integration, 1 or 2-page exam forms. Exam findings trigger codes, plans, outbound documents. Robust query tool. Free text type and touchscreen friendly. Two-way eRx sends/imports prescribed medications. Cloud or local server. Unmatched customer support <30 seconds, 98% satisfaction. Image management module works with any EHR.	Global Vision Equipment Booth 3171 772-281-2200 www.globalvisionequipment.com
503-707-8600 www.first-insight.com	GlobeChekBooth 318772-299-1404 www.globecheck.com
Fischer Surgical, Inc. Booth 440 314-303-7753 www.fischersurgical.com	BEST OF SHOW VIDEOS Check out the <i>AAO 2018 News</i> Sunday edition for information on this year's winners. View the videos at the Videos on Demand com- puter terminals in South Hall A and the Lakeside Center.



AMERICAN ACADEMY OF OPHTHALMOLOGY®

Visit the Academy Resource Center

Booth 508

Talk to experts about MACRA and MIPS

Demo the Academy's Ophthalmic News and Education (ONE®) Network and IRIS® Registry (Intelligent Research in Sight)

Get expert practice management and coding advice

Receive late-breaking federal and state governmental affairs and advocacy information

Renew your Academy, ISRS or AAOE membership – or join

ATTEND SPECIAL EVENTS:

Academy Foundation Donor Reception

Saturday, Oct. 27, 3 – 5 p.m. Museum of Vision, Booth 704 As a special thank-you, the Academy Foundation invites its donors to meet Academy leaders, learn about the future home of the Museum of Vision and enjoy refreshments.

*Ophthalmology** *Glaucoma* and *Ophthalmology** *Retina* Launch Celebration

Sunday, Oct. 28, 10 - 11 a.m. Resource Center, Booth 508 Subscribers, reviewers and published authors of the new Ophthalmology Glaucoma and Ophthalmology Retina journals are invited to meet editors-in-chief Henry D. Jampel, MD, and Andrew P. Schachat, MD, and the editorial boards.

EyeCare America[®] Volunteer Reception

Sunday, Oct. 28, 3 – 4:30 p.m. Museum of Vision, Booth 704

EyeCare America volunteers, join your colleagues as we honor you and your dedication to this vital public service. Residents who have pledged to volunteer on entering practice are also invited to learn firsthand how satisfying the experience can be.

2019 Focal Points® Spotlight Social

Monday, Oct. 29, 4 – 5 p.m. Resource Center, Booth 508

Subscribers, editors, authors and contributors to *Focal Points* are invited to preview and celebrate the new format, which launches January 2019. Each new *Focal Points* issue will examine current research and discuss practical ways to implement it. Find Resources to Enhance Your Clinical and Business Success



Save 10% on purchases*

Restrictions apply.

Good-Lite Company Booth 545 Hawaiian Moon Booth 662 727-709-1494 www.aloecream.biz **GOOD-LITE®** www.good-lite.com Heidelberg Engineering Booth 2236 508-530-7900 www.heidelbergengineering.com HEINE Booth 321 Good-Lite offers more than 2,000 products for acuity, low vision, color vision, ocular motility, contrast sensitivity, disease 800-367-4872 www.heine.com detection, binocular vision, pediatrics and more! We are the most comprehensive solution to your ophthalmic handheld testing and evaluation needs. Hablamos español! Stop by booth 545 and ask about our show specials! Highlights of Ophthalmology Booth 3120 507-301-0496 www.jphmedical.com 800-362-3860 www.good-lite.com GreenSky LLC Booth 3276 Howard Instruments, Inc. Booth 2612 844-810-7713 www.GreenSky.com 205-758-9083 www.howardinstruments.com Booth 3176 **Guangdong Fortune Newvision** Huai'An Chaoyue Medical Supplies Factory Booth 1971 Optoelectronics, Ltd 86-0159-6156-1416 www.chaoyue-medical.com 86-0769-3901-2155 fortune-newvision.com Booth 1268, 4362 **Guardion Health Sciences** Booth 1153 Huaian Frimen Co., Ltd 800-873-5141 www.guardionhealth.com 86-5178-989-0096 www.frimen.com Gulden Ophthalmics Booth 523 Hurricane Medical Booth 3403 215-884-8105 www.guldenophthalmics.com 941-751-0588 www.hurricanemedical.com Haag-Streit Group Booth 2226 Huvitz Co., Ltd Booth 4514 82-10-2486-7402 www.huvitz.com 513-398-3937 www.haag-streit-usa.com Hager Bambach Booth 239 HVO, Inc. Booth 1680 828-456-4455 www.hvoinc.com 828-624-3310 www.hagerbambach.com Booth 4057 HAI Laboratories, Inc. i.com medical GmbH Booth 2778 781-862-9884 www.hailabs.com 49-89520-328-80 www.icom-medical.de Hangzhou Mocular Medical Technology, Inc. Booth 153 IBM Watson Health Booth 767 86-0571-8187-5029 www.mocular.net 262-367-0700 www.merge.com Hausted Booth 3571 Icare USA, Inc. Booth 4600 770-368-4700 www.grahamfield.com 888-422-7313 www.icare-usa.com

IDxBooth 2977319-248-5620 www.eyediagnosis.net	International Council of Ophthalmology (ICO) Booth 1100 415-521-1651 www.icoph.org
Illuminix Industries Booth 2474 800-290-5140 www.illuminixeye.com Industries	IRIDEXBooth 3721650-940-4700 www.iridex.com
Ilooda Co., Ltd Booth 2048 82-22-6294-557	IrisVisionBooth 2172855-207-6665 https://irisvision.com
I-MED Pharma Booth 2409 I-514-685-8118 www.imedpharma.com	Ivantis, Inc.Booth 2273949-600-9650 www.ivantisinc.com
ImprimisRxBooth 4254858-704-4040 www.imprimisrx.com	Jaggi Ophthalmic Instrument Co.Booth 452191-1145-0526-692 www.jaggijaggi.com
Inami & Co., Ltd Booth 2800 81-338-14-5916 www.inami.co.jp/en/	JAMA NetworkBooth 559312-464-5000 www.jamanetwork.com
Indo Webal Surgical Pvt LtdBooth 65191-3324-3389-97 www.indowebal.com	Jaypee Highlights Medical Publishers, Inc.Booth 3120507 301-0496 www.jphmedical.com
Indo-German Surgical CorporationBooth 73291-2224-2218-09 www.indogerman.com	JCAHPO Booth 1505 800-284-3937 www.jcahpo.org
Infinity Massage ChairsBooth 537603-347-6006 www.infinitymassagechairs.com	JEDMED Booth 2200 314-845-3770 www.jedmed.com
Informational PostersBooth 1371415-561-8500 www.aao.org	Jenoptik Optical Systems, LLCBooth 4462561-881-7400 www.jenoptik.com
Innova Systems, Inc.Booth 2147630-920-8880 www.innovasystemsusa.com	JM SurgicalsBooth 241, 276586-2152-8188-36 www.jmsurgicals.com
Integra LifeSciencesBooth 4262901-417-7868 www.integralife.com	Johnson & Johnson VisionBooth 3708714-247-8200 www.jjvision.com
Integrated Orbital Implants, Inc.Booth 3119858-677-9990 www.ioi.com	Kai Industries Co., LtdBooth 204781-575-28-5040 www.kai-ind.com/global/en

Kala Pharmaceuticals, Inc.Booth 2576, 2676781-996-5252 www.kalarx.com	Labtician Ophthalmics, Inc.Booth 6211-905-829-0055 www.labtician.com
Karger PublishersBooth 2171860-675-7834 www.karger.com	Lacrimedics, Inc.Booth 2758253-964-0360 www.lacrimedics.com
Katalyst Surgical, LLCBooth 2662636-536-5950 www.katalystsurgical.com	LacriScience, LLCBooth 450301-523-4338 www.lacripen.com
KatenaBooth 1020973-989-1600 www.katena.comKeeler Instruments, Inc.Booth 2244	Lacrivera Booth 2613 LACRIVERA, a division of Stephens Instruments, is proud to offer the NEW VERAPLUG™ FLEX- FIT™ Punctal Occluder for the
610-353-4350 www.keelerusa.com KeraLink International Booth 2356 410-752-3800 www.keralink.org	treatment of dry eye. The FlexFit nose technology is designed to collapse when inserted allowing for easy sizing, unparalleled reten- tion and excellent patient comfort. Also available is the VeraPlug Punctal Occluder, Vera90 Ex- tended Wear Plug and VeraC7 Collagen Visit Poeth #2017 for
Khosla Surgical Industries Booth 2000 91-2223-8752-68 www.khoslasurgicals.com	Collagen. Visit Booth #2613 for more information. 855-857-0518 www.lacrivera.com
King Khaled Eye Specialist HospitalBooth 521966-11-482-1234 EXT 4749 www.kkesh.med.sa	Learning Lounge - American AcademyBooth 126of Ophthalmology415-447-0365 www.aao.org
Kirwan Surgical Products LLCBooth 2500781-834-9500 www.ksp.comImage: Comparison of the second sec	Leica MicrosystemsBooth 3703800-248-0123 www.leica-microsystems.com
Konan Medical USABooth 2515949-576-2100 www.konanmedical.com	Leiters Booth 1705 800-292-6772 www.leiterrx.com
Koryo Eyetech Co., LtdBooth 43982-25-48-1881 www.eyelikelens.com	LENSARBooth 1539888-536-7271 www.lensar.com
Kowa American CorporationBooth 3000800-966-5692 www.kowa-usa.com	Lenstec, Inc.Booth 2857727-571-2272 www.lenstec.com
Kugler Publications Booth 2821 31-2068-4-5700 www.kuglerpublications.com	Lifeline Medical Devices Pvt LtdBooth 195091-2402-6221-44 www.lifelineiol.com

LIGHTMED Booth 3468 949-218-9555 www.lightmed.com
LightScalpelBooth 558866-589-2722 www.lightscalpel.com
Lions Clubs International FoundationBooth 1777630-468-6791 http://www.lcif.org/EN/index.php
Lions Eye Institute For Transplant & Research Booth 413 813-289-1200 www.lionseyeinstitute.org
Lite Source, Inc. Booth 2846
502-266-9645 www.litesourceinc.com
Live Eyewear, Inc. Booth 553
Manufacturer of polarized fitover sunglasses, low vision aids and post-op eyewear. Featured brands include COCOONS, COCOONS LOW VISION, and LIGHTGUARD. See us for exclusive show specials!
800-834-2563 www.cocoonseyewear.com
LKC Technologies, Inc. Booth 3204
301-840-1992 www.lkc.com
Lombart Instrument Booth 2208
/5/-855-8888 www.iompartinstrument.com
Lumenis Vision Booth 3239 408-764-3000 www.lumenis.com
408-764-3000 www.iumenis.com
Luneau Technology USA Booth 3268 800-729-1959 www.luneautechusa.com
M&S Technologies Booth 2949
847-763-0500 www.mstech-eyes.com
Macan Medical Booth 2276 302-645-8068 www.macanmedical.com

This is NOT a plug.

This is a solution

to the discomfort of dry eye. The **VeraPlug**[®] **FlexFit**,[®] is simple to size and insert, comfortable for your patients, and comes with the confidence of our industry leading **45 day retention pledge**.

Visit us at AAO booth #2613

VeraPlug™ **FlexFit**™

lacrivera.com (855) 857-0518 2500 Sandersville Rd · Lexington KY 40511 USA

A FRESH PERSPECTIVE™

MacuLogix

MANI

Booth 2049

717-914-4067 www.maculogix.com

Madhu Instruments Pvt Ltd Booth 1849

91-011-4270-1030 www.madhuinstruments.com

Mallinckrodt Pharmaceuticals

908-238-6600 www.mallinckrodt.com

Booth 4203

Booth 1905

Booth 1810

Booth 1273

314-849-7773 www.crestpointmgt.com

Marasco & Associates, Healthcare Architects & Consultants

Marasco & Associates is a national architecture/consulting firm. We are the industry leaders developing 375+ ASC's and 1,750+ healthcare facilities nationwide including 300+

303-832-2887 www.mahca.com

Ophthalmic Facilities.

Marco Ophthalmic

TRS-5100 Refraction System

- Automated Digital Refractor • Cuts refraction time by
- 25-50%
- Replaces a standard refractor
- Programmable refractions
- Instantaneous comparison of old vs. new Rx
- Integrates data to other automated devices through Marco Connect and EMR
- TRS-5100 Product/Model name: **REFRACTOR RT-5100**

904-642-9330 www.marco.com

EYEPLAY EXPERIENCE

Take a break from AAO 2018 courses at the Academy's EyePlay Experience, available Saturday through Monday at Booth 2581. Check it out to:

- Give back to the community.
- Recharge your mobile device.
- Relieve stress with therapy animals.
- Challenge a colleague to Ping-Pong.
- Enjoy a complimentary seated massage.
- Take a selfie with a giant #AAO2018 hashtag.
- Get assistance at the Tech Bar.

Marco Ophthalmic

OPD-SCAN III

Corneal Topographer, Integrated Wavefront Aberrometer, Pupilometer and Autorefractor/Keratometer

- 23 diagnostic tests in 10 seconds/eye
- Separates corneal & internal astigmatism
- Corneal spherical aberrations and coma
- Retro illumination images
- Angle kappa and Angle Alpha
- White to white corneal size
- Toric IOL Summary marks axis pre-op • ECCP - Effective Central Corneal Power for IOL calculations

Integrated wavefront refractions with digital refractors & lens meters

904-642-9330 www.marco.com

Marco Ophthalmic

Booth 1810

ULTRA M SERIES Slit Lamps Series

- Light-gathering Ultra Optics, enhanced by high-luminance LED
- Comfortable 6 degree converging tubes
- Compact & vertical tower light sources
- Patent-pending, built-in background illumination
- Up to five magnification settings
- Blue-correction LED, red-free, yellow, and cobalt filters
- Larger 0-14 mm, continuously variable aperture
- Joystick with illumination boost or image capture trigger
- · Integrated, cable-free, low-voltage optimized for digital imaging integration
- Expanded chin rest for easier patient access

904-642-9330 www.marco.com

Marco Ophthalmic

iON IMAGING[™] Anterior Segment Imaging

- Intra-optics design
- Dedicated slit lamp imaging App
- On-board display for still and video images
- All controls at your fingertips
- Diagnose, image, integrate & educate in seconds
- Powerful educational teaching tool for staff and patients

904-642-9330 www.marco.com

CHECK OUT EYENET CORPORATE EVENTS

Leave room in your schedule for EyeNet's complimentary corporate educational events. Visit aao.org/eyenet/corporate-events for timing and other details.









Marco Ophthalmic

Boo<u>th 1810</u>

Booth 1810

EPIC® Refraction System

- Sets the standard for time efficient refractions in only a 4'x6' footprint
- Five technologies at a single position workstation
- Customized with 4 Objective Refraction systems (or Slit Lamp)
- Plus acuity system and choice of lensmeters
- Motorized adjustable table
- Color touch screen that can be turned to patient for near testing
- Perfect for delegation or simply most efficient refractive productivity

904-642-9330 www.marco.com

Marco Ophthalmic

TS-310 Tabletop Refractive Workstation

- Combined Digital Refractor and Chart System
- Ambidextrous Controls and dial button with S/C/A mode key
- Intuitive, high precision color LCD touch screen
- Same charts for far and near vision measurement
- Contrast sensitivity testing at 3 levels
- Rapid, quiet and smooth lens changes
- Flexible standing or seated measurement position
- · Versatile, compact footprint with simple installation

904-642-9330 www.marco.com

Market Scope Booth 2257 314-835-0600 www.market-scope.com

917-830-7960 www.marsdenadvisors.com

605-341-4595 www.mastel.com

Marsden Advisors

Mastel Precision

Meccanottica Mazza S.r.L Booth 2822 39-3629-0218-5 www.meccanotticamazza.it

Med X Change, Inc. Booth 2673 941-794-9977 www.mexchange.com

86-22-8371-3808 www.meda.com.cn

MedDev Corporation

800-543-2789 www.meddev-corp.com

MedEdicus

MEDA Co., Ltd

203-493-8077 www.mededicus.com

Medennium, Inc.

949-789-5340 www.medennium.com

Medical Ministry International Canada 1-905-664-0804 www.mmicanada.ca

Medical Technical Products

949-551-4762 www.medtechprods.com

Medicel AG

41-71-727-1050 www.medicel.com

Medinstrus UAB

370-685-86206 www.medinstrus.lt

Med-Logics, Inc.

949-582-3891 www.mlogics.com

941-359-3129 www.medone.com

MedOp Health, Inc.

MedOne Surgical, Inc.

Booth 453

888-290-6294 EXT 111 www.maxivision.com

Meds In Motion Pharmacy

Merge, an IBM Company

262-367-0700 www.merge.com

Booth 2578

Booth 767

801-506-6999 www.medsinmotion.com

2018 EXHIBITOR GUIDE • 21



Booth 2378

Booth 3019

Booth 1577

Booth 3068

Booth 934

Booth 562

Booth 3475

Booth 2139

Booth 865

Booth 2003

Booth 2404

Booth 3505



EyeNet Corporate **Lunches**

EyeNet® Magazine helps you make the most of your time at AAO 2018 by bringing you free corporate educational program lunches* onsite at McCormick Place.

Room E353c, Lakeside McCormick Place **Check-in and Lunch Pickup** 12:15-12:30 p.m. Lunches are provided on a first-come basis. **Program** 12:30-1:30 p.m.

Programs

Saturday, Oct. 27 Diabetic Eye Disease: Clinical Challenges and Practical Tips for Multidisciplinary Disease Management

Speakers: Harold Bays, MD (endocrinologist), John W. Kitchens, MD

Presented by Regeneron Pharmaceuticals, and designed for U.S. retina specialists.

Sunday, Oct. 28 INSiiGHTS AT AAO: A Spotlight on Dry Eye Treatment

Speakers: Eric D. Donnenfeld, MD, Edward J. Holland, MD, Terry Kim, MD

Presented by Shire

Monday, Oct. 29 Cataract Surgery: Life is Beautiful When the Pupil Behaves

Speakers: Eric D. Donnenfeld, MD, Cynthia A. Matossian, MD, FACS, Steven M. Silverstein, MD, Denise M. Visco, MD,


Keith A. Walter, MD Presented by Omeros Corporation, and designed for U.S. cataract surgeons.

Check aao.org/eyenet/corporate-events for updated program information.

* These programs are non-CME and are developed independently by industry. They are not affiliated with the official program of AAO 2018 or Subspecialty Day. By attending a lunch, you may be subject to reporting under the Physician Payment Sunshine Act.

MetrovisionBooth 184433-3-201-7195-0 www.metrovision.com	Museum of Vision (Foundation of the Booth 704 American Academy of Ophthalmology)
Micro Medical Devices, Inc. Booth 962 818-222-3310 www.micromedinc.com	AMERICAN ACADEMY OF OPHTHALMOLOGY* Protecting Sight. Empowering Lives. Visit the Museum of Vision to see "Remedy," an exhibit that features books on ophthalmic medicine published during the
MicroSurgical TechnologyBooth 4400425-556-0544 www.microsurgical.com	Renaissance and early modern periods. The 2018 exhibit, "Remedy," is devoted to the symptoms and treatments described in the Museum of Vision's rare book library. Display cases explore early descriptions of disease and
MID Labs, Inc. Booth 2849 510-357-3952 www.midlabs.com	physicians' use of natural plants, minerals, and metals as medi- cation from the 16th to early 20th centuries. 415-447-0315 www.faao.org
MIRA, Inc. Booth 3404 508-278-7877 www.mirainc.net	MynosysBooth 753510-770-9296 www.mynosys.com
Miracles In Sight Booth 1871 336-765-0932 www.miraclesinsight.org	Myvision Medical Pvt LtdBooth 286291-9909-6198-61 www.myvisionmedical.com
Modernizing MedicineBooth 1349561-880-2998 www.modmed.com	Nadia International, Inc.Booth 2608800-237-1868 www.ronadro.com
Modulight, Inc. Booth 162 358-207-43900-0 www.modulight.com	Natus Medical IncorporatedBooth 4248800-303-0306 www.natus.com
MORIA Booth 3733 33-1-467-4467-4 www.moria-surgical.com	Neitz Instruments Co., LtdBooth 150381-03-3237-0552 www.neitz-ophthalmic.com
MSI Precision Instruments Booth 2962 610-935-8181 www.msiprecision.com	NeoMedix CorporationBooth 1303714-258-8355 EXT 115 www.neomedix.net
MTI, Inc. Booth 1165 888-520-4998 www.mti.net	Netgain Technology, LLCBooth 2277877-797-4700 www.netgaincloud.com
WATCH AAO 2018 AT HOME Want free livestreams of approximately 20 hours of educational content from AAO 2018? Videos are available for viewing until Jan.	Network Technology Development, LtdBooth 1671852 6225-6941 www.networkmedical.com.hk
31, 2019. Access the Virtual Meeting with your Academy login and password at aao.org/virtual-meeting. Want to own your favorite presentations? Use AAO Meetings on	New Vision, Inc. Booth 237

Want to own your favorite presentations? Use AAO Meetings on Demand to watch full-motion videos, including audio and screen capture, from AAO 2018, Subspecialty Day, and the AAOE Program. 886-22-593-2028 www.nvi-laser.com



OD-OS GmbH

949-415-7255 www.od-os.com

Oertli Instrumente AG

41 71 747 4200

41-71-747-4200 www.oertli-instruments.com

Ohvem Sales Corporation

Booth 353

Booth 3003

91-9300-9400-95 www.aaa.org.in

Omeros Corporation

Booth 3772

Booth 3053

206-676-5000 www.omeros.com

OMIC-Ophthalmic Mutual Insurance Company Booth 1004

OMIC OPHTHALMIC MUTUAL INSURANCE COMPANY A Risk Retention Group

Visit the OMIC Exhibit and see our new sample patient education materials, informed consent documents, and other loss prevention tools.

OMIC is the largest insurer of ophthalmologists in the United States and is endorsed by more than 50 ophthalmic state and subspecialty societies, including the exclusive sponsorship of the American Academy of Ophthalmology.

800-562-6642 www.omic.com

Omni Lens Pvt Ltd Booth 1703

91-7927-5440-35 www.omnilens.in

Oncospark, Inc.	Booth 2774
601-329-0076 www.onospark.com	

OPHTEC

31-5052-5-1944 www.ophtec.com

ACADEMY RESOURCE CENTER

See the latest products and learn what services the Academy has to offer at the Resource Center (Booth 508).

Academy staff members are on hand at the Information desk and throughout the exhibit to answer your questions and help you zero in on the resources that will be most useful for your practice. In addition, *AAO 2018 News* will be distributed in the literature bins.



probes for endophotocoagulation since 2001. Our Bending Laser Probes, with or without widefield illumination, allow the surgeon to adjust and fix to an angle direction between 0 and 90 degrees. We can connect to all operating room lasers and light sources.

Call us for a no cost evaluation at 1 877 4 25 GAUGE (4284)

1-877-4-25-GAUGE (4284) www.ophthalmed.com

Ophthalmic Instrument Company

800-272-2070 www.oic2020.com

Ophthalmic Instruments, Inc.

Booth 2656

Booth 1951

800-373-1924 www.oiica.com

Ophthalmic Photographers' Society

Booth 2556



OPS is a professional organization of more than 1,000 members worldwide who share an interest in photography and imaging technology as applied to ophthalmology. Founded in 1969, the OPS remains true to its mission of providing its members an effective means to learn and advance the science and art of ophthalmic imaging.

417-725-0181 www.opsweb.org

Ophthalmology	Dispensing	Services	Boo
---------------	------------	----------	-----

201-787-0986 www.opticaldispensing.com

Ophthalmology Management

Booth 1808

Booth 2504


oth 338

215-628-6550 www.ophthalmologymanagement.com

Ophthalmology Times

732-346-3067 www.ophthalmologytimes.com





Optelec Booth 4414, 4 800-444-4443 www.vfogroup.com	415 Paragon BioTeck, Inc. Booth 734 866-424-1192 www.paragonbioteck.com
Optimetrics, Inc. Booth 32 305-863-4144 www.optimetrics.com	235 Pelion Surgical Booth 2905 803-649-9284 www.pelionsurgical.com
OPTIUSA/BK Frames Booth 23 305-685-0003 www.optiusa.com	Booth 122 41-41-552-4700 www.peschketrade.com
Optopol Technology Sp. z o.o. Booth 7 326-709-173 www.optopol.com	I34 PHAKOS Booth 1504 33-1-498-8740-0 www.phakos.com 33-1-498-8740-0 www.phakos.com
Optos, Inc. Booth 32 800-854-3039 www.optos.com	262 Phoenix Technology Group Booth 4244 925-485-1100 www.phoenix-clinical.com
Optovue, Inc. Booth 10 510-623-8868 www.optovue.com	D26 Physician Search Associates Booth 2973 919-656-4117
Orbis International Booth 20 646-674-5500 www.orbis.org	045 PhysIOL Booth 3419 32-43-61-0549 www.physiol.eu
Oyster Point Pharma Booth 12 609-382-9032 www.oysterpoinrx.com	278 Pine Pharmaceuticals Booth 3071 844-218-4138 www.pinepharmaceuticals.com
P.S.T. Corp. Booth 26 413-447-8051 www.pst-corp.com	629 PinPoint Optics Booth 2757 203-874-0448 www.pinpointoptics.net
Pacific Ophthalmic Supply Booth 18 650-637-9391 www.pacific-ophthalmic.com	B50 Plus Surgicals Booth 1171 91-8820-3059-93 www.plussurgicals.com
Pan-American Association Booth 12 of Ophthalmology (PAAO) 817-275-7553 www.paao.org	200 Plusoptix, Inc. Booth 1875 800-488-6436 www.plusoptix.com
Panasonic Industrial Medical Imaging Booth 1 201-392-6833 www.panasonic.com/visionsystems	169 PopUp Restaurant Experience - PURE Booth 2381 415-561-8500 www.aao.org
NEW EXHIBITOR PAVILION	Practice Flow Solutions Booth 832

678-935-7911 www.practiceflowsolutions.com

NEW EXHIBITOR PAVILION

The New Exhibitor Pavilion features all the latest industry offerings at AAO 2018 in 1 place. Visit Booths 2270-2879.

Precision VisionBooth 2412815-223-2022 www.precision-vision.com	Rectangle HealthBooth 2976914-298-2072 www.rectanglehealth.com
Premier Ophthalmic Services, Inc.Booth 1938800-597-7152 www.premierop.com	Refocus GroupBooth 563214-368-0200 www.refocus-group.com
Prescott's, Inc. Booth 2456 719-481-3353 www.surgicalmicroscopes.com	Regeneron Pharmaceuticals, Inc.Booth 107914-847-7000 www.regeneron.com
Pricon-Iscon Surgicals, Ltd Booth 2001 91-2912-7405-53 www.pricon.co.in	Reichert, Inc.Booth 3726716-686-4500 www.reichert.com
PRN Physician Recommended Nutriceuticals Booth 737 800-900-2303 www.prnomegahealth.com	Rendia Booth 2726 877-321-5481 www.rendia.com
Promedica International Booth 2358 714-460-7363 www.promedica-intl.com	RET, Inc. Booth 2467 82-10-3004-8672 www.retinc.kr Image: Second Secon
Qioptiq, An Excelitas Technologies Company Booth 3062 855-382-2677 www.excelitas.com	Review of OphthalmologyBooth 3162610-492-1000 www.reviewofophthalmology.com
Quantel MedicalBooth 3056877-782-6835 www.quantel-medical.com	RGP, Inc. Booth 1840 800-522-9695 www.rgpergo.com
Quark Pharmaceuticals, Inc.Booth 4162510-402-4020 www.quarkpharma.com	RICA Surgical ProductsBooth 3271847-671-5600 www.ricasurgical.com
Quest Medical, Inc.Booth 2945972-390-9800 www.questmedical.com	Rose Micro SolutionsBooth 522, 650, 1841, 2559716-608-0009 www.rosemicrosolutions.com
Quidel Corporation Booth 3073 Quidel's Eye Health solutions, InflammaDry and AdenoPlus, are rapid, lateral-flow based, point-of-care products for the	Royal BeeBooth 2674702-209-0462
detection of infectious and inflammatory diseases of the eye. 800-874-1517 www.quidel.com	Rumex InternationalBooth 2826727-535-9600 www.rumex.net
Ray Vision, Inc.Booth 2964949-656-0516 www.ray-vision.com	S4OPTIKBooth 2939888-224-6012 www.s4optik.com

Sanford HealthBooth 2473701-234-6510 www.sanfordhealth.org	Shanghai MediWorksBooth 2847Precision Instruments Co., Ltd86-21-5426-0421 www.mediworks.com.cn
Santen, Inc.Booth 4120415-268-9100 www.santenusa.com	Shanghai Xiaolai Science & Trade CompanyBooth 66486-2154-5400-28 www.rengongjingti.com
Scarf KingBooth 2071646-232-6105	Shenyang Xingqi Pharmaceutical Co., LtdBooth 65786-4006-1355-99 www.sinqi.com
SCHWIND eye-tech-solutions GmbHBooth 221549-6027-5-080 www.eye-tech-solutions.com	Shenzhen MOPTIM Imaging Technique Co., Ltd Booth 2667 86-7558-408-4505 www.moptim.com
SciCanBooth 1162800-572-1211 www.scican.com	ShireBooth 144617-349-0200 www.shire.com
ScienceBased HealthBooth 4054888-433-4726 www.sciencebasedhealth.com	Sight SciencesBooth 4317, 4711877-266-1144 www.sightsciences.com
Scrubs TailoredBooth 2575213-559-7566 www.scrubstailored.com	Sihong Ophsurin Co., LtdBooth 86686-5278-82-1082 www.chinaophsurin.com
sea-change, inc.Booth 2609561-629-3707 www.sea-change.com	SLACK Booth 2921 856-848-1000 www.healio.com/ophthalmology
Second Sight Medical ProductsBooth 540818-833-5060 www.secondsight.com	SMR Ophthalmic Pvt LtdBooth 204091-2228-7593-04 www.smrindia.com
SEE InternationalBooth 1838805-963-3303 www.seeintl.org	Soderberg InsrumentsBooth 2626800-222-8095 www.walmaninstruments.com
Segal Optiks Pvt LtdBooth 35291-9820-0277-45 www.segal.co.in	SolutionreachBooth 441866-826-5914 www.solutionreach.com
Senju Pharmaceutical Co., LtdBooth 107881-6-6201-9570 www.senju.co.jp	Sometech, Inc.Booth 416382-2-2025-1000 www.sometech.com
Shanghai Canton Optics Equipment Co., Ltd Booth 3076 86-2158-5889-77 www.chinacanton.com	Sonogage, Inc.Booth 4403216-464-1119 www.sonogage.com



Surgical Specialties Corp.

877-991-1110 www.surgicalspecialties.com

Booth 564

SRS Health

201-802-1300 www.srs-health.com

Booth 2845

SurgilumBooth 1946910-202-2202 www.surgilum.com	The Body Doctor, LtdBooth 126944-14848-6881-6 www.the-body-doctor.com
Surgistar, Inc.Booth 3854760-598-2480 www.surgistar.com	The Eye GroupBooth 2100561-852-9998 www.theeyegroup.com
Suzhou Kangjie Medical, Inc.Booth 316986-1526-2488-992 www.szkjmedical.com	The Lebanon CorporationBooth 2408800-428-2310 www.honanballoon.com
Suzhou Microclear Medical Instruments Co., Ltd Booth 3471 86-5126-7067-163 www.microcleartech.com	The OphthalmologistBooth 33944-15657-4-5200 www.theophthalmologist.com
Takagi Seiko Co., LtdBooth 113181-269-22-4512 www.takagi-j.com	THEA Booth 2532 33-4-739-8143-6 www.laboratoires-thea.fr
Tarun EnterprisesBooth 247691-9335-1545-56 www.optitecheyecare.com	Thieme Medical PublishersBooth 2904800-782-3488 www.thieme.com
Taylor & FrancisBooth 2173215-625-8900 www.tandfonline.com	ThromboGenics Booth 4063 32-16-751-310 www.thrombogenics.com
TBI / Tissue BanksBooth 2356410-752-3800 www.keralink.org	Tianjin Suowei Electronic Technology Co., Ltd Booth 4233 86-2283-7127-45 www.suowei.com.cn
TEACBooth 2872323-726-0303 http://teacmv.jp/en	Titan Medical, LLC Booth 1771 7-903-307-8572 www.titanmedical.su
Tear Film InnovationsBooth 2974844-458-9776 www.tearfilm.com	Titan SurgicalBooth 4205310-293-5757 www.titansurgical.com
Tecfen MedicalBooth 2050805-967-1153 www.tecfenmedical.com	Tomey CorporationBooth 273981-525-81-5327 www.tomey.com
Technology Pavilion - American Academy Booth 168 of Ophthalmology 415-447-0365 www.aao.org/2018	Topcon Medical SystemsBooth 3745201-599-5150 www.topconmedical.com
Tekia, Inc.Booth 1601949-699-1300 www.tekia.com	Tracey TechnologiesBooth 2508281-445-1666 www.traceytechnologies.com

TruForm Optics

800-792-1095 www.tfoptics.com

Booth 313

Booth 116

Booth 314

Booth 3219

TTI Medical 925-553-7828 www.ttimedical.com

630-572-2288 www.ubs.com/team/tategroup

UFSK-International OSYS, GmbH

UBS Financial Services

49-94178-862-15 www.ufsk-osys.com

United Sight	Booth 534
877-439-3226 www.unitedsight.org	
University of Wisconsin Department of Ophthalmology	Booth 2876
608-262-1334 www.eyes.wisc.edu	
US IOL, Inc.	Booth 3717, 4117
000 751 7010	

800-354-7848 www.usiol.com

USRetina	Booth 2864
844-698-7738 www.usretina.com	
Varitronics, Inc.	Booth 2145

610-356-3995 www.varitronics.com

Vector Vision	Booth 1153
800-873-5141 www.guardionhealth.com	
Verana Health	Booth 953
650-223-5560 www.digisight.net	
Virtual Field	Booth 2478

760-613-4364 virtualfield.io

Vision Associates Booth 1804 732-564-1100 www.visassoc.com














Buy Early • Save Money Limited Inventory Booth 545

Vision Group HoldingsBooth 4522561-965-9110 www.lasikvisioninstitute.com	Wilson OphthalmicBooth 931800-955-6544 www.wilsonophthalmic.com
Vision ShareBooth 341888-657-4448 www.visionshare.org	WinFame USA, Inc.Booth 550626-442-8238 www.winfameusa.com
VisionCare, Inc.Booth 4264408-872-9393 www.centrasight.com	Wolters Kluwer Booth 3300 215-521-8300 www.lww.com
Vital Tears Booth 3503 859-259-4925	World Glaucoma AssociationBooth 302031-20-575-4239 www.worldglaucoma.org
Vmax Vision Booth 1847 321-972-1823 www.vmaxvision.com	Zabby's Booth 141 91-11-2341-3805 www.zabbys.net
Volk Optical, Inc.Booth 2244440-942-6161 www.volk.com	ZeaVision LLC Booth 1948 314-628-1000 www.eyepromise.com
Vortex SurgicalBooth 3064636-778-4350 www.vortexsurgical.com	Zeiss Booth 3736 925-557-4100 www.zeiss.com/us
VRmagic, Inc. Booth 1544 617-444-8761 www.vrmagic.com	Ziemer Ophthalmics Booth 1056 618-462-9301 www.ziemergroup.com
Walman Instruments/ Eye Care AllianceBooth 2626800-222-8095 www.walmaninstruments.com	Zocular Booth 4062 844-965-8527 www.zocular.com
Welch Allyn Booth 312 315-685-4100 www.welchallyn.com	Zomex Instruments Co. Pvt LtdBooth 106892-52-353-0266 www.zomexinstruments.com
West TeleVox SolutionsBooth 2653800-644-4266 www.televox.com	Zoomax (USA), Inc. Booth 2777 866-296-8388 www.zoomaxusa.com
Wicab, Inc. Booth 152 608-829-4500 www.wicab.com	Zumax Medical Co., LtdBooth 337686-512-6690-3235 www.zumaxmedical.com
Wills Eye HospitalBooth 113215-928-3000 www.willseye.org	NEED A BREAK? VISIT A LOUNGE. In Grand Concourse, Level 3: YO, SO, and International Lounges. In Room S104b: AAOE Member Lounge.



Floor plan is accurate as of Sept. 7, 2018.

Learn about Interventional Glaucoma **at AAO 2018**

Ellex technology, which includes advanced SLT Selective Light Therapy and ABiC[™] canalbased glaucoma surgery, performed with the iTrack[™] surgical system, gives you the ability to intervene across the full spectrum of the glaucoma disease process – and enables significant numbers of patients to minimize or obviate the need for anti-glaucoma medications. Join us at AAO 2018 Exhibit #1044 for a series of physician-led initiatives to see first-hand how Ellex's interventional glaucoma technologies can help you to prioritize proactive glaucoma care.

MEET THE EXPERT	MYTH BUSTER
SAT, 27 OCT	SAT, 27 OCT
SUN, 28 OCT	SUN, 28 OCT
10:00 - 11:00am	10:30 - 11:00am
3:00 - 4:00pm	3:00 - 3:30pm



MAHMOUD KHAIMI, MD



Meds

ABiC



LEARN MORE AT WWW.ELLEX.COM/AAO

*These sessions are not affiliated with the official program of AAO 2018

