

QCDR Measure:

IRIS-6: Acquired Involutional Entropion: Normalized lid position after surgical repair

National Quality Strategy Domain:

Effective Clinical Care

Measure Type:

Outcome

Description:

Percentage of surgical entropion patients with a postoperative normalized lid position.

Instructions:

This measure is to be reported **each time** a entropion surgical procedure is performed during the reporting period. The measure is intended to reflect the quality of services provided by the surgeon performing the procedure.

Denominator:

All patients aged 18 years or older with a diagnosis of involutional entropion who underwent a surgical procedure for the condition.

Denominator Criteria**Patients aged ≥ 18 years**

AND

Diagnosis of Involutional Entropion

ICD-9 [for use 1/1/2015 – 9/30/2015]

- Senile entropion (ICD-9: 374.01)

ICD-10 [for use 10/1/2015 – 12/31/2015]

- Senile entropion (ICD-10: H02.031, H02.032, H02.033, H02.034, H02.035, H02.036, H02.039)

AND

Entropion Surgical Procedure

- Lateral tarsal strip (CPT: 67917)

- Lateral canthoplasty (CPT: 67950)
- Entropion repair (CPT: 67924)
- Excision and repair of eyelid (CPT: 67966)
- Excision and repair of eyelid (CPT: 67961)

Numerator:

Patients with normalized lid position postoperatively

Numerator Options:

Performance Met:	Patients who achieved normalized lid position postoperatively within 90 days following surgery
Performance Not Met:	Patients who did not achieve normalized lid position postoperatively within 90 days following surgery

Improvement Notation:

Higher score indicates better performance